



45807C/Revised: February 2008

MURI-LUBE[®] *LIGHT MINERAL OIL*

DESCRIPTION:

Muri-Lube is a mixture of liquid hydrocarbons obtained from petroleum, and is composed entirely of saturated aliphatic and naphthenic hydrocarbons. It is stabilized by the addition of Vitamin E. Muri-Lube is completely free from color, odor and taste, and is non-drying, non-staining, not gummy or sticky. It is sterile and non-supportive to pathogenic bacteria. It has excellent lubricating, penetrating and softening properties.

INDICATIONS AND USAGE:

As a lubricant for surgical instruments. Muri-Lube is not intended for parenteral or oral use.

HOW SUPPLIED:

Product Number	Order Number	Vial Size	Mineral Oil
25402	6332325402	2 mL	Light
25410	6332325410	10 mL	Light

Packaged in 25s.

Store at 20° to 25°C (68° to 77°F) [see USP Controlled Room Temperature].

Vial stoppers do not contain natural rubber latex.



APP Pharmaceuticals, LLC
Schaumburg, IL 60173

45807C
Revised: February 2008