

 $4\,5\,8\,0\,7\,D$  /Revised April 2025

## **MURI-LUBE**® LIGHT MINERAL OIL

## DESCRIPTION:

**DESCRIPTION:**Muri-Lube is a mixture of liquid hydrocarbons obtained from petroleum, and is composed entirely of saturated aliphatic and naphthenic hydrocarbons. It is stabilized by the addition of Vitamin E. Muri-Lube is completely free from color, odor and taste, and is non-drying, non-staining, not gummy or sticky. It is sterile and non-supportive to pathogenic bacteria. It has excellent lubricating, penetrating and softening properties. properties.

INDICATIONS AND USAGE:
As a lubricant for surgical instruments.
Muri-Lube is not intended for parenteral or oral use.

HOW SUPPLIED:
Product Order
Number Number
25402 632325402
632325410 Mineral Oil Light Light Vial Size 2 mL 10 mL Packaged in 25s.

Store at 20° to 25°C (68° to 77°F) [see USP Controlled Room Temperature].

The container closure is not made with natural rubber latex.

www.fresenius-kabi.com/us

45807D

Revised: April 2025